

Power Chairs Reference no. PC013 Cruiser 12A



The Cruiser 12A has been introduced specifically to overcome restrictions which relate to motorised transport on P&O and Cunard cruise ships. Where clients are unable to obtain accessible cabins to take a mobility scooter the cruiser 12A can be hired. The chair folds like a manual wheelchair for easy storage. On board are 2x 11.6Ah lithium ion batteries which can give a max speed of 6.4kmh(4mph) The range and speed of the chair may vary depending on the weight of the user, terrain type, outside temperature and battery charge. We advise that you plan your route to avoid adverse terrain, inclines, kerbs and steps if possible; if this is unavoidable, we do not advise an incline above 5° and a kerb / step above 4cm (1.5"). If it has been agreed that adjustments are necessary, this will need to be done prior to 1pm on the day of the cruise. Please arrange this with the office.

Max user weight	100kg / 15.7st	Kerb climber	No
Total weight	45kg	Seat raising	No
Total width	64cm	Seat recline	No
Total height	96cm	Backrest recline	No
Seat width	45cm	Footplate / standard leg rest / elevating leg rest	Standard leg rest
Seat height	49cm	Footrest length	TBC
Seat depth	43cm	Seat angle	N/A
Backrest height	47cm	Backrest angle	N/A
Armrest height	21cm	Calf supports	No (Available on request)
Turning diameter / radius	ТВС	Flip up armrest	No
Ground clearance	6cm	Headrest	No
Foldable / comes apart	Yes / (removable leg rests)	Swivel seat	No
Removable batteries	No	Drive wheel position	Rear wheel
Anchor points / crash tested	No	Battery Amp hour	8.8Ah (Version 1) 11.6Ah (Version 2)
Control options	User joystick (L/H or R/H)	Battery type	2x 24V Lithium
Swing-away joystick	No	Tyre type / colour	Solid / Grey
Extra manual brakes	Yes / No (please specify if required)	Max Speed	4mph
Lights and hazards fitted	No	Length	100 cm

Disclaimer