

# MobilityatSea

## Power Chairs

Reference no. PC009

Cruiser 9 - Mirage Power Chair



The Mirage power chair is suited to the less dependent user. The padded seat and armrests provide comfort whilst the stability of rear wheel drive makes this one of the easier chairs to manoeuvre. On board are 2x 35Ah AGM dry cell batteries which can give a max speed of 6kmh(3.8mph) and a max range of 19km(12miles), although the range and speed of the chair may vary depending on the weight of the user, terrain type, outside temperature and battery charge. We advise that you plan your route to avoid adverse terrain, inclines, kerbs and steps if possible; if this is unavoidable we do not advise an incline above 5° and a kerb / step above 4cm(1.5"). If it has been agreed that adjustments are necessary, this will need to be done prior to 1pm on the day of the cruise. Please arrange this with the office.

## Specifications

Ground Clearance:	9.5cm / 3.75"
Length:	107cm / 42.5"
Width:	61cm / 24.5"
Seat Width:	45.75cm / 18"
Turning Radius:	57cm / 22.5"
Weight with batteries:	72kg / 158.5lbs
Maximum user weight:	120kg / 264lbs

## Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.