



Wheelchairs

Reference no. WC019

Bariatric Transit Wheelchair



This wheelchair is designed for the larger passenger. The brakes are positioned on the handles for attendant control

We advise that you plan your route to avoid adverse terrain, inclines, kerbs and steps (if this is unavoidable, we do not advise an incline above 5 degrees and kerb or step above 3cm/1.1").

Please be mindful that gangways and airbridges can vary in the degree of gradient due to the tides, some products are not designed to cope with acute slopes.

THIS PRODUCT NEEDS TO BE FOLDED AND PUSHED THROUGH NARROW CABIN DOORWAYS.

Specifications

Height	101cm/40"
Width	73.7cm/29"
Closed Width	35.6cm/14"
Length	85.1cm/33.5"
Seat Height	50.8cm/20"
Seat Width	60cm / 24"
Seat Depth	46cm / 18"
Weight	14.6kg / 32.3lbs / 2.3st
Max User Weight	178kg / 392lbs / 28st

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which it is intended.