



Scan for 'how to' video



Mobility Scooters

Reference no. ST020

Stowaway Colibri



The Stowaway Colibri disassembles into 5 pieces for easy storage. The range and speed of the scooter may vary depending on the weight of the user, terrain type, outside temperature and battery charge. We advise that you plan your route to avoid adverse terrain and inclines if possible; if this is unavoidable, we do not advise an incline above 5°.

Please be mindful that gangways and air bridges can vary in the degree of gradient due to the tides, some products are not designed to cope with acute slopes.

Specifications

Length	101cm / 40"
Width	50.5cm / 20"
Height	Up to 80.5cm / 31.6"
Total Weight	49kg / 109lbs / 7st7
Heaviest Piece	15.6kg / 35lbs / 2st4
Max User Weight	120.6kg / 266lbs / 19st
Turning Radius	220cm / 87"
Battery Type	SLA - dry cell
Removable Battery	Yes

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.