



Scan for 'how
to' video



Powerchairs

Reference No. PC041

Dashi MG



The Dashi MG provides all the comfort and manoeuvrability of an electric wheelchair. Please be advised that the range and speed of the chair may vary depending on the weight of the user, type of terrain, outside temperature and battery charge. We advise that you plan your route to avoid adverse terrain, inclines, kerbs and steps (if this is unavoidable, we do not advise an incline above 5 degrees and kerb or step above 3cm (1.1")).

Please be mindful that gangways and air bridges can vary in the degree of gradient due to the tides, some products are not designed to cope with acute slopes.

THIS PRODUCT NEEDS TO BE FOLDED AND PUSHED THROUGH NARROW CABIN DOORWAYS.

his chair performs best on the 3-bar speed setting, using the joystick to control the speed. We do not advise the slowest setting (1-bar) unless you are a confident user as this setting may impair performance.

Specifications

Length	106cm / 41.7"
Width	60cm / 23.6"
Height	91.44cm / 36"
Folded Width	34cm / 13.3"
Seat Width	45cm / 18"
Seat Depth	43.18cm / 17"
Weight	17.6kg / 38.8lbs / 2st7
Max User Weight	101kg / 222lbs / 16st
Battery Type	Lithium Ion

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.