



Moving and Handling

Reference no. MH016

Reliant 350 Stand Aid



The Reliant 350 can be used for standing and seated transfers. It has manual adjustable legs that splay open for closer transfers. The user must be able to weight bear. It is the client's responsibility to provide their own loop attachment sling.

Specifications

Weight (product)	36kg (79.5lbs / 5.6st)
Width (between arms)	(max) 74cm (29 inches)
Width (overall)	93cm (36.6 inches)

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which it is intended.