



Moving & Handling

Reference No. MH022

Alerta MoveAssist Transfer Aid



This Transfer Aid has manual adjustable legs that splay open for closer transfers. **The user must be able to stand up and remain standing until the paddles are lowered into position and vice versa.**

Specifications

Overall Length	94cm / 37"
Overall Height	110cm / 43.3"
Height of Seat (min)	67.8cm / 26.6"
Height of Seat (max)	87.6cm / 34.5"
Height of Base Plate	10cm / 3.9"
Height of Knee Pad (min)	29.5cm / 11.6"
Height of Knee Pad (max)	56cm / 22"
Width Between Legs (min)	60cm / 23.6"
Width Between Legs (max)	79.5cm / 31.2"
Seat Width	51cm / 20"
Overall Width (min)	66cm / 26"
Overall Width (max)	92cm / 36.2"
Weight	36kg / 78.4lbs / 5.6st
Max User Weight	190.5kg / 420lbs / 30st

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.