



Moving & Handling

Reference No. MH016

Reliant 350 Stand Aid



The Reliant 350 can be used for standing and seated transfers. It has manual adjustable legs that splay open for closer transfers. **The user must be able to weight bear.**

It is the client's responsibility to provide their own loop attachment sling

Specifications

Overall Length	109cm / 43"
Overall Height (max)	168.9cm / 66.5"
Height of Base Plate	7.6cm / 3"
Height of Knee Pad (max from base)	58.4cm / 23"
Height of Spreader Bar (max)	161.5cm / 63.5"
Leg Width (min)	55cm / 21.6"
Leg Width (max)	84.5cm / 33"
Width between arms (max)	74cm / 29"
Overall Width	93cm / 36.6"
Weight	36kg / 79.5lbs / 5.6st
Max User Weight	152.4kg / 336lbs / 24st

Disclaimer

Certain circumstances may determine that another product is supplied. We will endeavour to notify you in such an event however this will not impact on the functionality of the product or the purpose of which is it intended.